# Stem Cells Generation



# 5-Azacytidine

DNA hypomethylation reagent that increases stem cell reprogramming efficiency.<sup>1</sup> Enhances the efficiency of iPS cell generation.<sup>2</sup>

<u>Product No: 10-1267</u> 50 mg 250 mg

**BAY K8644** 

Enables reprogramming of Act3/Klf4-transduced embryonic fibroblasts in combination with BIX-01294.<sup>3</sup>

Product No: 10-4518 5 mg 25 mg

Cordycepin

AMPK activator that maintains stem cell pluripotency and increases iPS cell generation efficiency.<sup>4</sup>

Product No: 10-2408 10 mg 50 mg

**DZNep** 

Methyltransferase inhibitor that can be used in chemical cocktails to induce iPSs from somatic cells.<sup>5</sup>

Product No: 10-1524 1 mg 5 mg

Na Valproate

Improves stem cell reprogramming efficiency and enables efficient induction of pluripotency without introduction of the oncogene c-Myc.<sup>6</sup>

<u>Product No: 10-1009</u> 5 g

PD 325901

Induces generation of ground state iPS cells from human and rat somatic cells in combination with CHIR-99021 and A83-01.7

Product No: 10-2161 5 mg 25 mg

Pifithrin α

p53 inhibitor that induced pluripotency in primary bone marrow derived mesenchymal stromal cells in conjunction with SB431542 and PD325901.8

Product No: 10-2480 5 mg 25 mg

**Purmorphamine** 

Part of a chemical cocktail that reprograms mouse somatic cells into pluripotent stem-like cells<sup>9</sup> and human astroglial cells into functional neurons<sup>10</sup>.

Product No: 10-4622 5 mg 25 mg

RepSox

Potent ALK5 inhibitor that can replace Sox2 in reprogramming adult cells into pluripotent stem cells.<sup>11</sup>

Product No: 10-4546 5 mg 25 mg

**Thiazovivin** 


Dramatically improves the efficiency of induced pluripotent stem cell generation from human fibroblasts.<sup>12</sup>

Product No: 10-1191 5 mg 25 mg

Tranylcypromine hemisulfate

Induces reprogramming of human primary keratinocytes into iPS cells in conjunction with CHIR99021.<sup>13</sup>

Product No: 10-4535 50 mg 250 mg



# **Purmorphamine**

### Thiazovivin

### REFERENCES

- 1. Mikkelsen et al. (2008) Nature 454 49
- 2. De Carvalho et al. (2010) Trends Cell Biol. 20 609
- 3. Shi et al. (2008) Cell Stem Cell 3 568
- 4. Wang et al. (2020) Sci. Rep. 10 2187
- 5. Hou et al. (2013) Science 341 651
- 6. Hangfu et al. (2008) Nat. Biotechnol. 26 795
- 7. Li *et al.* (2009) Cell Stem Cell **4** 16
- 8. Megges et al. (2015) Stem Cell Res. 15 565
- 9. Kang et al. (2015) Biomaterials 35 7336
- 10. Zhang et al. (2015) Cell Stem Cell 17 735
- 11. Ichida et al. (2009) Cell Stem Cell 5 491
- 12. Lin et al. (2009) Nat. Methods 6 805
- 13. Li et al. (2009) Stem Cells 27 2992
- 14. Wang et al. (2019) Cell Biol. 97 702

400 Davis Dr. Suite 600 Plymouth Meeting, PA 19462 610-994-1134 sales@focusbiomolecules.com focusbiomolecules.com