Ubiquitin-Proteasome

Proteasome Inhibitors

MG-132

Reversible inhibitor of the 20S proteasome (K_i = 4 nM).¹ Most commonly used proteasome inhibitor. Stimulates autophagy.

Product No: 10-1309

5 mg 25 mg_

Bortezomib (Velcade)

Potent, selective and reversible proteasome inhibitor ($K_i=0.6$ nM). ² Inhibits proliferation of a number of tumor cell lines ($IC_{50}=7$ nM).

Product No: 10-2120 5 mg 25 mg

Lactacystin

A Streptomyces metabolite that is thought to bind irreversibly to the active site N-terminal threonine residue of the catalytic β -subunit of the 20S proteasome, thereby inhibiting its chymotrypsin- and trypsin-like activities.³

Product No: 10-1314

200 µg 1 mg

Epoxomicin

Epoxomicin is a potent, selective and irreversible inhibitor of the 20S proteasome.⁴ It does not inhibit papain, chymotrypsin, trypsin, calpain or cathepsin B at 50 μ M.⁴

Product No: 10-1384

Carfilzomib (Kyprolis)

An analog of epoxomicin, Carfilzomib is a selective proteasome inhibitor currently approved as a therapeutic for multiple myeloma.⁵

100 µg

Product No: 10-1477 5 mg 25 mg

Deubiquitinase (DUB) Inhibitors

P5091

Selective inhibitor of USP7 (IC_{50}=4.2 μM). 6 Induces apoptosis in multiple myeloma cells.

Product No: 10-1422 5 mg 25 mg

PR-619

A nonselective, reversible inhibitor of DUBs (5-20 $\mu M)^7$ that is highly useful for preserving ubiquitinated proteins during cell lysis.

Product No: 10-1485 5 mg 25 mg

OTUDin3

Potent inhibitor of the deubiquitinase OTUD3, a key enzyme in tumorigenesis, neurodegeneration, metabolism, and immunity.⁸

Product No: 10-4981 5 mg 25 mg

iBAP-II

Selective Inhibitor of the H2AK119 deubiquitinase BAP1.9

Product No: 10-4038 5 mg 25 mg

WP-1130

A novel DUB inhibitor of USP9x, USP5, USP14 and UCH37.¹⁰ WP1130 induces degradation of ERG, an important driver of prostate oncogenesis overexpressed in 40% of prostate tumors.¹¹

Product No: 10-1351 5 mg 25 mg

Check out our extensive collection of Molecular Glues!



DUB and Proteasome Inhibitors

Use PR-619 to preserve ubiquitinated proteins in cell lysates

DUBs become unregulated in cell lysates and will act quickly to deubiquitinate your protein of interest during cell lysis. To preserve maximum amounts of ubiquitinated protein, include PR-619 (50-100 μ M) in your lysis buffer.⁷

REFERENCES

- 1. Tsubuki et al. (1996), J. Biochem., **119** 572
- 2. Adams et al. (1999), Cancer Res., 59 2615
- 3. Fenteany et al. (1995), Science, 268 726
- Meng *et al.* (1999), Proc. Natl. Acad. Sci. USA, **96** 10403
- 5. Kuhn et al. (2007) Blood, **110** 3281
- 6. Chauhan et al. (2012) Cancer Cell, 22 345
- 7. Altun et.al. (2011) Chemistry & Biology,18 1401
- 8. Zhang et al. (2023) Cell Death Dis., 14 378
- 9. Tsuboyama et al. (2022) Oncogene 41 2152
- 10. Kapuria et al. (2010) Cancer Res., 70 9265
- 11. Wang et al. (2014) Proc Natl Acad Sci U S A., 111 4251

400 Davis Dr. Suite 600 Plymouth Meeting, PA 19462 610-994-1134 sales@focusbiomolecules.com focusbiomolecules.com