


Catalog # 10-2204 Cisplatin

CAS# 15663-27-1
Cis-Diaminodichloroplatium
CDDP; cic-Diamineplatinum(II) dichloride
Lot # X106257



Induces intra- and interstrand DNA adducts which poison topoisomerase I.¹ Synergizes with topoisomerase I inhibitors such as topotecan.¹ Induces apoptosis by sequential activation of caspase-8, -3 and -6 in osteosarcoma cells.² Inhibits expression of X-linked inhibitor of apoptosis protein (XIAP) in human LNCaP cells.³ Numerous mechanisms inhibit propagation of cisplatin-induced DNA damage to signal apoptosis resulting in resistance.⁴ Clinically useful anticancer agent.

- 1) Van Waardenburg et al. (2004), Platinated DNA adducts enhance poisoning of DNA topoisomerase I by camptothecin; J. Biol. Chem,, **279** 54502
- 2) Siddik et al. (2003), Cisplatin: mode of cytotoxic action and molecular basis of resistance; Oncogene, 22 7265
- 3) Seki et al. (2000), Cisplatin (CDDP) specifically induces apoptosis via sequential activation of caspase-8, -3 and -6 in osteosarcoma; Cancer Chemother. Pharmacol., **45** 199
- 4) Nomura et al. (2004), Cisplatin inhibits the expression of X-linked inhibitor of apoptosis protein in human LNCaP cells; Urol. Oncol., **22** 453
- 5) Raghavan et al. (2015), Dimethylsulfoxide inactivates the anticancer effect of cisplatin against myelogenous leukemia cell lines in in vitro assays.; Indian J. Phamracol., **47** 322

PHYSICAL DATA

Molecular Weight: 300.05

Molecular Formula: Cl₂H₆N₂Pt

Purity: 98% by HPLC

Elemental Analysis: (Conforms)

Solubility: Soluble in Water (up to 2.5 mg/ml with warming)

Note: Compound will dissolve in DMSO, but DMSO inactivates cytotoxicity.⁵

Physical Description: Yellow solid

Storage and Stability: Store as supplied, desiccated at room temperature for up to 1 year from the date of purchase.

Solutions in distilled water are not stable and must be used within one working day. Do not store

solutions.

Materials provided by Focus Biomolecules are for laboratory research use only and are not intended for human or veterinary applications.