




Catalog # 10-4548

PAPA NONOate

146672-58-4

N-[3-Aminopropyl(propyl)amino]-N-hydroxynitrous amide

Lot # FBA1281



PAPA NONOate is a nitric oxide donor and dissociates to the free amine and NO in a pH dependent manner. The half-life of PAPA NONOate in 0.1 M phosphate buffer at pH 7.4 is 15 minutes at 37 °C, or 77 minutes at room temperature.¹ Dissociation is essentially instantaneous at a pH of 5.0.

- 1) Keefer *et al.* (1996), "NONOates" (1-substituted diazen-1-ium-1,2-diolates) as nitric oxide donors: Convenient nitric oxide dosage forms; *Methods Enzymol.*, **268** 281

PHYSICAL DATA

Molecular Weight:	176.22
Molecular Formula:	C ₆ H ₁₆ N ₄ O ₂
Purity:	97% by (¹ H-NMR)
Solubility:	Water or Ethanol
Physical Description:	Off-white solid
Storage and Stability:	Store as supplied at -80°C for up to 1 years from the date of purchase. Protect from exposure to air or moisture. Solutions are not stable and must be prepared immediately prior to use. See half life information above.

Materials provided by Focus Biomolecules are for laboratory research use only and are not intended for human or veterinary applications.

Focus Biomolecules LLC 400 Davis Drive, Suite 600 Plymouth Meeting PA 19462

www.focusbiomolecules.com